
High-end benefits in entry-level digital







Experience a world of high-end benefits with Neo






What features do you need to solve your client's problems?

Until now, it was always necessary to use high-end digital hearing instruments to get the features that you want for your clients.

Neo, Connected to the world

Renew the connection between your client and their world.

- Open fitting with large vents
- □ Adaptive feedback cancellation
- Noise reduction
- Directional microphones
- Multiple programs
- Multiple channels
- **Expansion**
- Automatic telecoil
- Telecoil boost
- Push button control
- **FM compatibility**
- **Reliable and accurate software**





High-end features



High-end features are now accessible to everybody.

Neo has all of the essential elements to bring satisfaction to your client. Features that you once could only have in high-end hearing instruments are now available in the entry-level category.

Neo has:

- ✓ OpenFit[™]
- ✓ Adaptive Feedback Cancellation
- ✓ Adaptive Noise Reduction in 8 frequency bands
- **Directional microphones with active** microphone matching
- Multiple listening programs for different environments
- 🗹 5 channels for maximum flexibility
- ✓ Soft Noise Management[™]
- **M** Beep indicators for important functions
- ✓ Automatic telecoil
- **V** Telecoil boost
- **V** Push button control of hearing instrument functions
- Volume Control
- **F**M compatibility








Open up a new world of fitting possibilities with the perfect combination of Neo and OASIS plus.

The high performance fitting software platform to ensure a perfect fit every time.



Intelligent software makes your job easier.

- **V** First fit in as few as 3 mouse clicks
- ✓ 5 channel fitting accuracy with 5 MPO controls
- **V** Fine-tuning controls at your fingertips
- ✓ TriQualizer[™] for precise adjustment
- **V** 3-D view of hearing instrument response
- ✓ Interactive fine-tuning to support the dialogue with your client
- COSI for counselling and ensuring that your client's goals are met
- **MOAH** link compatibility
- ✓ Industry standard programming cables and FlexConnect
- ✓ Unified fitting flow with all Bernafon digital hearing instruments

Neo and OASIS plus, the perfect combination of advanced hearing instruments and high performance fitting software





Improved hearing in noisy situations

Directional microphones with active microphone matching, proven technology to enhance the perception of speech in noise

 Adaptive Noise Reduction in 8 frequency bands to improve listening comfort and reduce fatigue
Compatibility with ear level FM equipment such

as LEXIS

Optimal hearing for any environment

Multiple listening programs that can be easily customized for any need


 ✓ Soft Noise Management[™] to reduce the intrusiveness of every day environmental sounds
✓ NAL NL1 fitting rationale, adjusted to your client with the Client Experience Level



Comfort

✓ OpenFit[™] to reduce the occlusion effect and give a natural sensation of hearing

- Adaptive Feedback Cancellation, an automatic active feedback cancellation system to search for and eliminate annoying whistling
- ✓ 5 channels with independent controls to ensure that sound levels are comfortable



V European quality and design

- Easy to use program selection with a push button
- ✓ Intelligent beep indicators for hearing instrument functions








✓ Easy to fit and powerful fine tuning with the TriQualizer[™]

- Adaptive Feedback Cancellation without a reduction of high frequency gain, to preserve speech intelligibility
- ✔ OASIS plus has a proven track record of stability and reliability





Easy to explain




✓ Intelligence in the design and implementation of the features









Features that were once too expensive for many of your clients are now accessible.





A full range of hearing instruments for virtually any requirement

| Hearing Instrument | OSPL90, Peak / dB SPL [Earsim] | Full-On Gain, Peak/dB [Earsim] | Battery size | Local controls |
|-----------------------------|-----------------------------------|-----------------------------------|-----------------|----------------------------|
| Neo 102 BTE | 128 | 64 | 13 | Push button+Volume control |
| Neo 105 BTE DM | 127 | 66 | 13 | Push button+Volume control |
| Neo 112 BTE | 136 | 76 | 13 | Push button+Volume control |
| Neo 202 ITE | 129 | 62 | 13 | Push button+Volume control |
| Neo 352 HS (US/CAN only) | 124 | 56 | 312 | Push button+Volume control |
| Neo 322 ITC | 124 | 56 | 312 | Push button+Volume control |
| Neo 355 HS DM (US/CAN only) | 122 | 54 | 312 | Push button+Volume control |
| Neo 315 ITC DM | 122 | 54 | 312 | Push button+Volume control |
| Neo 302 ITC | 120 | 50 | 312 | Push button+Volume control |
| Neo 301 ITC | 120 | 50 | 312 | Push button |
| Neo 411 MC (Lab opt.) | 114 | 45 | 10 | Push button |
| Neo 401 CIC | 114 | 45 | 10 | Push button |







Bernafon, provides a wide range of digital hearing systems that cover all of your client's needs. All hearing instruments are fit with OASIS plus – the single, intelligent fitting software for all Bernafon instruments.



Innovative Hearing Solutions

www.bernafon.com

Bernafon AG Morgenstrasse 131 3018 Bern Switzerland Phone +41 (0)31 998 15 15 Fax +41 (0)31 998 15 90 www.bernafon.com

Bernafon Australia Pty. Ltd. Unit 12, 97 Castlemaine Street Milton QLD 4064 Australia Freecall 1800 809 111 Phone +61 (0)7 3250 0300 Fax +61 (0)7 3250 0372 www.bernafon.com.au Bernafon UK Cadzow Industrial Estate Low Waters Road, Hamilton ML3 7QE Scotland United Kingdom Phone 01 698 28 59 68 Fax 01 698 42 14 56 www.bernafon.co.uk

Bernafon New Zealand Ltd. Level 1, Building F 27-29 William Pickering Drive Albany, Auckland New Zealand Phone +64 9 920 42 40 Fax +64 9 920 42 43 www.bernafon.co.nz